dedicio®

CRP Test (test strip)

#311003D-10/311003D-20

DE	Gebrauchsanweisung	2
EN	Instructions for use.....	5
FR	Instructions d'utilisation.....	8
ES	Instrucciones de uso	11
IT	Istruzioni per l'uso	14
PL	Sposób użycia	17
CZ	Návod k použití.....	20
FI	Käyttöohje.....	23
SE	Användarinstruktioner	26
	Symbols	32

dedicio+
dedicated to life

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Spezifität: 98,7%

14. Referenzen

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417.
2. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
3. Macy EM, Hayes TE and Tracy RP (1997) Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. Clin. Chem. 43, 52-58.

Rev. 7, 2016-06-09 OM/TW

1. Intended Use

The dedicio® CRP Test is intended to be used for semiquantitative determination and monitoring of CRP concentrations in whole and capillary blood samples as well as in serum and plasma samples.

2. Introduction and Clinical Significance

C-reactive protein (CRP) in patient samples has been found in association with acute infections, necrotic conditions and a variety of inflammatory disorders. There is a strong correlation between levels of CRP and the onset of the inflammatory process. Monitoring the levels of CRP in patient samples indicates the effectiveness of treatment and the assessment of patient recovery. It is used in particular to differentiate bacterial infections from virus infections.

3. Principle of the Test

The dedicio® CRP Test is an immunochromatographic test based on two specific antibodies against human CRP. The concentration-dependent formation of test lines allows a rapid semiquantitative determination of CRP in whole blood/serum/plasma samples.

The tip of the test strip is dipped into the test sample which is diluted with buffer solution. The sample now migrates along the test strip from bottom to top. If the test sample contains CRP, it attaches to the first anti-CRP antibody which is conjugated to red colloidal gold for colour marking. The red CRP antibody-gold complex, together with the sample liquid, diffuses through the membrane that is pre-dispensed with lines of different amounts of the second anti-CRP antibody. The CRP-antibody-gold complex is immobilized by antibodies precoated on the membrane leading to the formation of red lines. The number of lines depends on the CRP concentration in the sample. The more CRP is contained in the sample, the more red lines become visible.

The control line serves as a procedural control and confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

4. Reagents and Materials Supplied

- 10/20 dedicio® CRP test strips (individually packaged)
- 10/20 plastic tubes containing buffer solution (1 mL)
- 10/20 end-to-end capillaries (10 µl)
- 1 workstation
- 1 package insert

5. Additional Required Materials

- Timer
- Lancet
- Capillary holder / tweezers

6. Storage & Stability

The test kit should be stored at 2-30°C until the expiry date printed on the sealed foil pouch. The test strip must remain in the sealed pouch until use. Do not freeze. Care should be taken to protect the components of the test kit from contamination. Do not use the test if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipments, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- This test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is, therefore, recommended that these products be treated as potentially infectious and handled observing usual safety precautions (e.g., do not ingest or inhale).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to testing.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Do not interchange or mix reagents from different lots.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

Preparation

Before performing the test, make sure that all test components are brought to room temperature (see "Storage & Stability"). Take a tube containing buffer solution out of the test kit. Mark it with patient name or ID. Open the screw cap.

Blood Sample Taking


1. Disinfect the finger tip. Using a lancet, extract a drop of blood from the finger tip.
2. With the supplied capillary tube, take a volume of 10 µl from the blood drop. It is important that the end-to-end capillary is filled until the upper end. Due to hygienic reasons, hold the capillary with a capillary holder or tweezers. Alternatively, the blood can also be taken with a micropipette.

Note: In case of using micropipettes or other capillaries a suitable volume of exactly 10 µl must be administered.

Dilute the blood sample immediately to avoid clotting.

Sample Dilution / Sample Stability



1. Administer the blood-filled end-to-end capillary into the corresponding plastic tube with dilution buffer. Alternatively, 10 µl of blood can be added directly with a micropipette into the buffer.

OR

transfer 5 µl of serum/plasma into the corresponding plastic tube with dilution buffer.

- Close the tube and shake the sample by hand forcefully for approximately 10 seconds so that the blood is released from the capillary and the sample and dilution buffer mix well.
- Let the diluted sample rest for approximately 1 minute.
- The sample can then be used immediately or stored for up to 8 hours.

Note: EDTA, citrate or heparin blood can be used as well. Before performing the test, it should be diluted accordingly with the supplied buffer.


9. Procedure of the Test

After opening the foil pouch the test should be performed quickly, as test strips are sensitive to humidity.

- Open the foil pouch and take out the test strip at the end labelled with CRP.

Do not touch the white result field with your fingers.

- Open the tube with the diluted sample and immerse the test strip with the other end into the liquid. Make sure to dip the test strip not above the maximum marking line (MAX). Avoid in any case a direct wetting of the result field by lateral dashes or by dipping the test strip to deep.
- Leave the test strip for a minimum of 10 seconds in the diluted sample until the slightly pink-coloured liquid front becomes visible in the result field.
- Take out the test strip and put it on a plain underlay that does not dehydrate the test strip (e.g. the test pouch). Alternatively, the test strip can remain in the tube. Start the timer.
- Wait for the appearance of coloured lines. Read the result after exactly 5 minutes. Please stick to this time limit to ensure correct semiquantitative results.

10. Interpretation of the Results

Please note that the test must be read after exactly 5 minutes.

To interpret the results please look at the line pattern that has formed in the white membrane field.


The cut-off of detection of the test is 10 µg/ml.

If more than one line (control line C) appears below the holder, the test result has to be regarded as positive: that CRP was detected in the sample.

Positive:


1 red test line (T1) and 1 red control line (C):

At CRP concentrations of approximately 10 µg/ml and higher, a test line T1 that is faint in the beginning appears in the lower area in addition to the control line.


2 red test lines (T1 and T2) and 1 red control line (C):

From approximately 40 µg/ml on a second faint coloured test line T2 appears above the first test line.


3 red test lines (T1, T2 and T3) and 1 red control line (C):

From approximately 80 µg/ml on a third faint test line T3 appears above the first two test lines.

Negative:

Only the red control line appears. The CRP concentration in the sample is <10 µg/ml.

Invalid:

The red control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for failure of the control line. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note: Please read the test only once after 5 minutes. It is quite natural for immunochromatography and the kinetics of rapid tests that the intensity of all lines increases over time. If the test is read later, an accurate semiquantitative interpretation of the test result is not ensured anymore, as low concentrations could then be estimated as too high.

11. Quality Control

The internal procedural control is included in the test strip. A coloured line appearing in the control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

12. Limitations

- CRP is not a specific marker for a certain disease. As with all *in-vitro* diagnostic devices, the result should not be interpreted on its own but in correlation with clinical findings. CRP increase often occurs before any symptoms become apparent, therefore, temporal connections should be considered as well.
- Intraindividual variation of CRP values is relatively high. In general, values >10 µg/ml can be regarded as elevated in the majority of the patients.
- The test is not suitable for risk stratification of coronary heart diseases. For this, it is not sensitive enough.

13. Performance Characteristics

The test was calibrated against the international human CRP reference standard of the WHO 85/506. The detection limit of the dedicio® CRP Test is adjusted to 10 µg/ml in such a way that a concentration of 10 µg/ml in the undiluted specimen material after dilution with the supplied buffer yields a positive red line.

If the undiluted specimen contains more than 40 µg/ml CRP, the second red line appears. At concentrations higher than 80 µg/ml CRP, the third red line appears.

These concentration ranges could be confirmed in comparison to a quantitative reference test.

Hook effects could not be proven up to a concentration of 2000 µg/ml CRP. At this concentration, the test indicated results of >80 µg/ml.

Specimen Correlation

The study was performed on 150 negative specimens (CRP<10 µg/ml, EIA confirmed, 150 plasma) and 233 positive serum and plasma specimens (CRP in the range of 10-40 µg/ml, EIA confirmed, 91 plasma; CRP in the range of 40-80 µg/ml, EIA confirmed, 60 plasma; and CRP >80 µg/ml, EIA confirmed, 82 plasma). All the results were determined at 5 minutes. The results are shown in the table:

Method		EIA			
		Negative	Positive		
dedicio® CRP Test	Ranges	0-10 µg/ml	10-40 µg/ml	40-80 µg/ml	>80 µg/ml
	0-10 µg/ml	148	2	0	0
	10-40 µg/ml	2	88	2	0
	40-80 µg/ml	0	1	57	4
	>80 µg/ml	0	0	1	78
Total results		150	91	60	82
% Agreement		98,7%	96,7%	95%	95,1%
			95,7%		

Sensitivity: 99,1%

Specificity: 98,7%

14. References

- Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417.
- Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
- Macy EM, Hayes TE and Tracy RP (1997) Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. Clin. Chem. 43, 52-58.

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Tarkkuus: 98,7%

14. Lähteet

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis IE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417.
2. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet:980-983.
3. Macy EM, Hayes TE and Tracy RP (1997) Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. Clin. Chem. 43, 52-58.

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nal von minden GmbH